

# A high-throughput platform to quantify bacterial defense against protist predation

Andrew M Hogan<sup>1,2</sup>, Alexandra Koumoutsi<sup>3</sup>, Nassos Typas<sup>1</sup>, Jordi van Gestel<sup>2</sup>

1. Molecular Systems Biology Unit, European Molecular Biology Laboratory, Heidelberg, Germany
2. Developmental Biology Unit, European Molecular Biology Laboratory, Heidelberg, Germany
3. Microbial Automation and Culturomics Core Facility, European Molecular Biology Laboratory, Heidelberg, Germany

## Background

- Protists are major predators of bacteria and thus a selection pressure towards the evolution of defense mechanisms, many of which are known<sup>1</sup> (Figure 1)
- "Microbial dark matter" (unannotated genes) may play roles in natural settings, such as interactions and predation
- No quantitative high-throughput assays exist to screen large strain libraries for predation defense
- Screening a large number of isolates may reveal new predation defense mechanisms that contribute to microbial ecosystem shaping





Figure 1. Prototypical bacterial defenses against protist predation

## Aims


1. Develop a fluorescence plate-based assay to screen for bacterial defense against *Dictyostelium discoideum* predation
2. Screen a global library of *Escherichia coli* isolates for predation defense
3. Validate the assay with gain and loss of function genetic screens

## Methods

### 1 Assay Design



### 2 Reporter Conjugation



### 3 Quantification




Figure 2. Protist predation assay design and quantification. "Tetrads" with duplicates of control and test prey strains are analyzed separately. The predator and prey are robotically pinned on SM/5 agar plates. *D. discoideum* spreads on the agar surface and consumes fluorescently labelled bacteria, the signals of which are recorded and processed in three ways to get rates of predation.

## Results

Total Colony Fluorescence Measurement Validly Identifies Resistant Isolates




Figure 3. Initial hits from screening *E. coli* library. A) Representative images from isolate plate 5 over time. B) The most resistant isolates are consumed up to 4-fold slower than the *E. coli* B/r control. C) Predation defense of *E. coli* ZG-35.1 shown in timelapse.

Different Predation Metrics Reveal the Diversity of Defense Phenotypes




Figure 4. Diverse defense strategies exist in the *E. coli* isolates. A) All vs all correlation suggests differences in how each metric scores defense. B) Some defense phenotypes are related to biofilm formation.

## Key Outcomes

- Developed and validated a high-throughput quantitative assay to screen for bacterial isolates with defense against protist predation
- A spectrum in strength and diversity of predation defense phenotypes was observed in a global collection of 553 *E. coli* isolates
- This assay will be useful for any bacteria-protist pairs that can grow on an agar surface

## Outlook

- Gain- and loss-of-function genetic screens to identify predation defense genes
- Examine effect of defense genes on single cell predator-prey interactions
- Screen libraries of soil bacterial isolates for resistant isolates

## Contact Information

andrew.hogan@embl.de  
@MicrobeAndrew  
vangestelgroup.com



We're curious to try other protist-bacteria combinations. If you have bacterial strain libraries or a model protist, and you'd like to find resistance mechanisms, reach out to us!

## References

1. Jousset A. 2012. Environ. Microbiol. 14(8): fiac057. DOI: 10.1093/femsec/fiac057
2. Schlechter et al. 2018. Front. Microbiol. 9: 3052 DOI: 10.3389/fmicb.2018.03052
3. Galardini et al. 2017. eLife. 6: e31035. DOI: 10.7554/eLife.31035

## Aknowledgements

Image Analysis:  
Christian Tischer and Arif Khan  
(EMBL Centre for Bio-Image Analysis)

